

CLINICAL MONITORING CHART FOR DENGUE FEVER

	Time & Date				
Parameter					
Sensorium	<ul style="list-style-type: none"> • clear & lucid • restless/delirious 	<ul style="list-style-type: none"> • clear & lucid • restless/delirious 	<ul style="list-style-type: none"> • clear & lucid • restless/delirious 	<ul style="list-style-type: none"> • clear & lucid • restless/delirious 	<ul style="list-style-type: none"> • clear & lucid • restless/delirious
Body temperature					
Pulse rate					
Blood pressure					
Pulse pressure (N) - > 20 mm					
Respiratory rate					
Pulse volume	<ul style="list-style-type: none"> • Good • Weak & thready • Feeble or absent 	<ul style="list-style-type: none"> • Good • Weak & thready • Feeble or absent 	<ul style="list-style-type: none"> • Good • Weak & thready • Feeble or absent 	<ul style="list-style-type: none"> • Good • Weak & thready • Feeble or absent 	<ul style="list-style-type: none"> • Good • Weak & thready • Feeble or absent
Capillary refill time	<ul style="list-style-type: none"> • ≤2s • >2s 	<ul style="list-style-type: none"> • ≤2s • >2s 	<ul style="list-style-type: none"> • ≤2s • >2s 	<ul style="list-style-type: none"> • ≤2s • >2s 	<ul style="list-style-type: none"> • ≤2s • >2s
Warmth of extremities	<ul style="list-style-type: none"> • Warm • Cold 	<ul style="list-style-type: none"> • Warm • Cold 	<ul style="list-style-type: none"> • Warm • Cold 	<ul style="list-style-type: none"> • Warm • Cold 	<ul style="list-style-type: none"> • Warm • Cold
Abdominal pain /vomiting					
Bleeding					
IVF infusion rate					
JVP					
Chest- Air entry, Basal creps					
Liver size					
Urine output					

Look for fluid overload as evidenced by

- Increase in heart rate >25 / mt.
- Increase in respiratory rate >5 / mt.
- Increase in liver size >2cm
- Filling up of JVP